# Instructions for Use

## TEMPORARY CB

Prepared: 13/07/2020 PRNT-106836 Rev 00

The following Instructions for Use are for Formlabs biocompatible photopolymer Temporary CB Resin (A2, A3, B1, C2 according to the VITA\* classical shade system). Basic information about safety and environmental concerns are also included. For more detailed safety and environment information, please refer to the Safety Data Sheet available at [dental.formlabs.com](https://www.formlabs.com). For further information regarding the use of this material, please contact Formlabs.

\*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

## 1. Intended purpose/Indication

Temporary CB resin is indicated for the fabrication of temporary dental restorations in conjunction with extra-oral light-curing equipment.

## 2. Contraindications

Known allergy to one or more ingredients. In cases of doubt, the allergy should be clarified and ruled out based on a specific test prior to the application of this device. Temporary CB Resin should not be used for purposes other than temporary crown and bridge work. Any deviation from these Instructions for Use can have negative effects on the chemical and physical quality of plastics manufactured from Temporary CB Resin.

## 3. Specific Manufacturing Considerations

### NOTIFICATION

The device specifications have been validated using the printer process parameters indicated below.

### REQUIREMENTS

All accessories used must be designated for Temporary CB Resin only. For full compliance and biocompatibility, Temporary CB Resin requires a dedicated resin tank, Stainless Steel Build Platform, Form Wash, and finishing kit which should not be mixed with any other resins.

## RECOMMENDED 3D PRINTER AND PRINTING PARAMETERS

- a. Hardware: Formlabs SLA 3D Printer
  - Laser wavelength: 405 nm
  - Form 2 Resin Tank LT or Form 3 Resin Tank
  - Stainless Steel Build Platform only
- b. Software: Formlabs PreForm
  - STL file import
  - Manual/Automatic rotation and placement
  - Manual/Automatic generation of supports
- c. Printing Parameters
  - Layer thickness: 50µm
  - Optimal orientation: horizontal orientation, occlusal plane facing the build platform
  - Crowns, inlays, onlays, veneers: Minimum wall thickness: 1.0mm
  - Bridges:
    1. Minimum wall thickness:
      - 1.0mm (anterior)
      - 1.5mm (posterior)
    2. • 12mm<sup>2</sup> (anterior)
      - 16mm<sup>2</sup> (posterior)
    3. Max. construction length for bridges 7 units
    4. Pontic max. one molar width
- d. Recommended Post-Processing Equipment:
  - Formlabs Form Wash
  - Isopropyl alcohol (IPA) ≥ 99%
  - Formlabs Form Cure
  - Sandblaster
  - Glass bead blasting material - Perlablast® micro\* ( BEGO Co.) recommended
  - Finishing kit
  - Dental composite polish and tools

\*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

## 4. Hazards And Precautions

### HAZARDS

Temporary CB Resin contains polymerizable monomers which may cause skin irritation, allergic skin reaction, serious eye irritation, respiratory irritation or other allergic reactions in susceptible persons. If resin contacts skin, wash thoroughly with soap and water.

If any discomfort or adverse health effects occur from contact with parts printed with Temporary CB Resin, discontinue use of this material.

Detailed information about the handling of Temporary CB Resin can be found in the Safety Data Sheet and [Formlabs.com](http://Formlabs.com).

1. **Eye Contact:** Can cause serious eye irritation. If eye irritation persists, seek medical attention.
2. **Skin Contact:** May cause skin irritation or allergic skin reactions. If skin irritation or rash occurs, seek medical attention.
3. **Inhalation:** May cause respiratory irritation. Avoid breathing mist, vapors, and sprays of this resin.
4. **Ingestion:** Avoid ingesting uncured or cured resin. If accidentally ingested, call a POISON CENTER and seek medical attention.
5. **Protection:** Protective clothing, glasses, and nitrile gloves should be worn while handling Temporary CB Resin.
6. **Environment:** May cause harmful effects to aquatic life. Avoid release into the environment.
7. **Disposal:** Dispose of contents and cartridge as per local regulations.

**Contains:** Esterification products of 4,4'-isopropylidiphenol, ethoxylated and 2-methylprop-2enoic acid, silanized dental glass, methyl benzoylformate, diphenyl(2,4,6-trimethylbenzoyl) phosphine oxide.

Total content of inorganic fillers (particle size 0.7 $\mu$ m) is 30-50% by mass.

### CAUTION

Wash the printed part with solvent in a properly ventilated environment with proper protective masks and gloves.

1. Expired or uncured Temporary CB Resin shall be disposed in accordance with local regulations.
2. IPA shall be disposed in accordance with local regulations.

## 5. Procedure for fabricating dental restorations with Temporary CB Resin

Shake well before using a new cartridge of Temporary CB Resin. Temporary CB Resin in the resin tank does not need to be mixed prior to printing, up to a month of non-use, as solids (fillers) in the material do not settle.

### A. PRINTING AND POST-PROCESSING

1. **Set up:** Insert resin cartridge, Stainless Steel Build Platform, and the compatible resin tank into the Formlabs 3D printer.
2. **Printing:**
  - a. Prepare a print job using Preform software. Import a dental restoration STL file. Orient files horizontally with the occlusal plane facing the build platform and generate supports. For recommendations on print orientation and support placement, see the detailed application guide at [support.formlabs.com](https://support.formlabs.com)
  - b. Send print job to printer. Begin print by selecting print job from print menu. Follow any prompts or dialogs shown on the printer screen. Printer will automatically complete the print.
3. **Washing:** Restorations printed with the Temporary CB Resin should be washed using the Form Wash as follows,
  - a. Remove the build platform with the printed parts from the printer and place it into a Form Wash with clean IPA ( $\geq 99\%$ ).
  - b. Set the Form Wash for 3 minutes. Make sure the parts are fully submerged in IPA when washing. Exceeding wash duration may affect dimensional accuracy and performance of printed parts over time.
  - c. Once parts are washed, inspect parts for any remaining uncured resin on print surfaces. Use an IPA-filled squeeze bottle to remove any remaining resin on the printed parts and in between supports and rafts. An IPA-soaked brush may be used to help with uncured resin removal.
4. **Drying:** Use compressed air to dry the parts. The dried parts will have a white, powdery coating on the print surfaces. If any wet, uncured resin is still present after drying, reuse the IPA-soaked brush and IPA-filled squeeze bottle to remove uncured resin and dry again.
5. **Part removal:** After drying, remove printed parts from the build platform by wedging the part removal tool under the print raft and rotating the tool. For detailed techniques, see the support material on [dental.formlabs.com](https://dental.formlabs.com).
6. **Post cure:** Restorations printed with Temporary CB Resin can be post cured in the Form Cure. To maintain dimensional accuracy and biocompatibility, post curing must be done in two steps:

- a. With the raft and supports still intact, place printed parts in the Form Cure with the raft side down.
- b. Form Cure parts at 60°C (140°F) for 20 minutes.
- c. Remove parts from Form Cure and use a handpiece with a cutting disc to separate supports and raft from printed parts.
- d. Carefully sandblast the printed part surfaces to remove the white, powdery coating using a sandblaster with glass bead blasting material 50 µm, such as Perlablast® micro\* (e.g. from BEGO Co. REF 46092/54302) sandblasting medium at a maximum blasting pressure of 1.5 bar.
- e. Check printed part fit and finish and contour parts as needed using a dental handpiece and carbide rotary burs.
- f. Place finished parts in the Form Cure and post cure again at 60°C (140°F) for an additional 20 minutes, with the occlusal plane facing upwards. For best results and full compliance, printed parts must be post cured twice as recommended. Any deviations from this protocol, may lead to unsatisfactory fit and affect biocompatibility.

**Attention:** Inspect the dental restorations for cracks. Discard if any cracks or damage are detected.

## B. POLISHING

Post cured parts must be polished using pumice stone and polishing compound prior to use.

1. High gloss polishing can be achieved with commonly used dental composite polishers.
2. Avoid overheating the parts during polishing.
3. Optimal surface quality is achieved by polishing after the full post curing cycle.
4. If rough marks are left behind after support removal, sand and polish to smooth surfaces to improve patient comfort.

\*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

## C. CLEANING AND DISINFECTION

Fully cured restorations made with Temporary CB Resin can be easily cleaned and disinfected with IPA or ethanol. Steam cleaning using e.g. Triton SLA\* (BEGO Co.), or MD 520\* (Dürr Dental Co.) impression disinfectants are also possible. Follow the manufacturer's instructions.

## D. LUTING AGENT COMPATIBILITY

The finished restoration can be attached to a patient's tooth using conventional temporary cements such as Temp Bond NE\* (Kerr Co.). If the subsequent permanent restoration will be attached with methacrylate-based composite cement, then eugenol-free temporary cements are recommended for use. Follow the Instructions For Use provided by the luting agent manufacturer.

## E. STORAGE

1. Store finished parts in closed, opaque or amber containers in a cool, dry place out of direct sunlight.
2. Store the resin cartridges at room temperature (approx. 22°C or 72°F). Do not exceed 28°C (82°F) or drop below 4°C (39°F) when in storage.
3. Avoid ignition sources.
4. Do not use resin cartridges past the expiration date on the cartridge label.

## F. DISPOSAL

Any cured resin is non-hazardous and may be disposed of in regular waste.

1. Follow clinic or facility protocols for waste that is considered as a biohazard.
2. Liquid resin should be disposed of in accordance with governmental regulations (community, regional, or national).
  - a. Contact a licensed professional waste disposal service to dispose of liquid resin.
  - b. As with all foreign substances, do not allow resin to enter in storm or sewer drainage systems.
  - c. Avoid release into the environment.
  - d. Contaminated packaging: Dispose of as unused product.

\*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

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Wilh. Herbst GmbH & Co. KG**
**Manufactured for: Formlabs Inc.**