INSTRUCTIONS FOR USE

Custom Tray Resin

EN	Instructions for Use	GA	Treoracha Úsáide
DE	Gebrauchsanweisung	HR	Upute za upotrebu
FR	Mode d'emploi	MT	Istruzzjonijiet għall-Użu
IT	Istruzioni per l'uso	LV	Lietošanas instrukcija
ES	Instrucciones de uso	LT	Naudojimo instrukcijos
PT	Instruções de uso	ET	Kasutamise juhised
NL	Gebruiksaanwijzing	SK	Návod na použitie
DA	Brugsanvisning	SL	Navodila za uporabo
SV	Bruksanvisning	RO	Instrucțiuni de utilizare
FI	Käyttöohjeet	PL	Instrukcja użytkowania
CS	Návod k použití	BG	Инструкции за употреба
HU	Használati útmutató	JP	使用説明書
GR	Οδηγίες χρήσης	ZH	使用说明

Prepared: 02/03/2023 PRNT-0016 Rev 02

INSTRUCTIONS FOR USE

Custom Tray Resin

Keep away from sunlight

Use-by date

Consult Instructions for Use

European Conformity

Batch Code

Catalog Number

Manufacturer

Temperature Limit

Authorized representative
in the European Community

Medical Device

Environmental hazard

Warning: Read Safety Information

Unique device identifier

Indicates the authorized
representative in Switzerland

Formlabs Ohio Inc.
27800 Lemoyne Rd,
Millbury, OH 43447
USA
+1 617 855 0762

EU Importer
Formlabs GmbH
Nalepastr. 18
12459 Berlin, Germany
+49 1573 5993322

Curaden AG
Amlehnstrasse 22
6010 Kriens
Switzerland

Before printing, refer to the Manufacturing Guide and dental applications guide for this material at support.formlabs.com for printing and post-processing recommendations and requirements to ensure the correct and safe usage of this material. For detailed safety and environmental information, the Safety Data Sheet is available at support.formlabs.com.

I. Introduction And Indications For Use

INTRODUCTION AND INTENDED PURPOSE

Custom Tray Resin is a light-curable polymer-based resin designed for the fabrication of biocompatible, short-term use, removable dental appliances by additive manufacturing. It is intended for 3D printing dental appliances such as dental impression trays.

INDICATIONS FOR USE

Custom Tray Resin is a photopolymer resin that can be utilized to print dental impression trays.

CONTRAINDICATIONS

Custom Tray Resin should not be used with a known allergy or hypersensitivity history to photopolymer resins. In case of an allergic reaction, stop use immediately, and consult a physician.

INTENDED PATIENT POPULATION

All patients and populations in which a clinician has prescribed a treatment associated with the intended purpose.

INTENDED USERS

Dental professionals who fabricate dental appliances.

CLINICAL BENEFIT

The clinical benefit of Custom Tray Resin is the ability to print dental appliances such as impression trays.

II. Safety

RESIN ASSOCIATED WARNINGS

Custom Tray Resin (uncured) contains polymerizable monomers.

- Please read Custom Tray Resin SDS at support.formlabs.com for hazards associated with liquid resin.

RESIN USE ASSOCIATED PRECAUTIONS

- Protective glasses and nitrile gloves should be worn while handling resin. For information about handling Custom Tray Resin please read Safety Data Sheet at support.formlabs.com.

INFORMATION FOR PATIENT SAFETY

- The patient should inform the clinician of any known allergy or hypersensitivity history to photopolymer resins.
- In case of an allergic reaction, stop use immediately and consult a medical physician.
- Do not use appliances if damage or cracks are present.
- If signs of extreme wear, deterioration, and/or perforations appear, the appliances should be replaced to avoid material failure and/or potential bodily injury.

III. Performance Characteristics

For detailed information on the mechanical properties of this material, technical data sheets can be found at support.formlabs.com. Custom Tray Resin has been evaluated in accordance with ISO 10993-1:2020 and ISO 7405:2018 and meets the requirements of all applicable biocompatibility testing.

IV. Specific Manufacturing Considerations

Refer to the Manufacturing Guide at support.formlabs.com for printing recommendations and requirements to ensure the correct and safe usage of this material.

- Hardware: Formlabs 3D Printer and accessories
- Software: Formlabs PreForm
- Part Orientation: Design guidelines for part orientation during printing can be found in the

Manufacturing Guide.

- Part Thickness: Design guidelines for part thickness can be found in Manufacturing Guide
- Recommended Post-Processing Equipment:
 - Formlabs validated accessories
 - Formlabs validated washing units with approved solvent
 - Formlabs validated curing units

REQUIREMENTS

For biocompatibility compliance, Custom Tray Resin was validated with a dedicated resin tank, build platform, wash unit and post-processing equipment that were not mixed with any other resins.

A. PRINTING

Printing guidelines and validated equipment can be found in the Manufacturing Guide at support.formlabs.com.

B. WASHING

Washing guidelines and validated equipment can be found in the Manufacturing Guide.

C. POST-CURING

Post-Curing guidelines and validated equipment can be found in the Manufacturing Guide.

D. SUPPORT REMOVAL

Support marks can cause abrasion if not removed and polished. See the Manufacturing Guide for polishing recommendations.

E. CLEANING & DISINFECTION

The appliances may be cleaned and disinfected according to facility protocols. Tested disinfection method: soaking the finished appliance in fresh 70% IPA for 5 minutes. Do not leave the part in the alcohol solution for longer than 5 minutes. After disinfection, inspect the part for damage or cracks to ensure that the integrity of the designed part meets performance requirements.

F. STORAGE

1. Printed appliances should be stored in closed, opaque, or amber containers. Excess light exposure over time may affect the color and performance of printed appliances.
2. Store cartridges at 10°C - 25°C (50°F - 77°F). Do not exceed 25°C (77°F) when in storage.
3. Keep cartridges closed & away from ignition sources.

G. DISPOSAL

1. Cured resin is non-hazardous and may be disposed of as regular waste.
2. Follow facility protocols for waste that may be considered biohazardous.
3. Unused liquid resin & solvents should be disposed of in accordance with local regulations.
4. Cartridge & contaminated packaging should be disposed of in accordance with local regulations.

H. ADVERSE EVENTS

To report any serious injury (temporary or serious permanent deterioration of health, a serious public health threat, or death), contact Formlabs Inc. (+1 617 855 0762) or your local Competent Authority.

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