# Instructions for Use

## PERMANENT CROWN

Prepared: 12/01/2021 PRNT-106965 REV 02

The following Instructions for Use are for Formlabs biocompatible photopolymer Permanent Crown Resin. Basic information about safety and environmental concerns are also included. For more detailed safety and environment information, please refer to the Safety Data Sheet available at [dental.formlabs.com](https://dental.formlabs.com). For further information regarding the use of this material, please contact Formlabs.

\*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

## 1. Indications For Use

Permanent Crown is indicated as an indirect restorative for both anterior and posterior restorations, including occlusal surfaces. The Permanent Crown material is used for fabricating permanent restorations such as inlays, onlays, veneers, and full crown restorations.

## 2. Contraindications

Known allergy to one or more ingredients. In cases of doubt, the allergy should be clarified and ruled out based on a specific test prior to the application of this product. Permanent Crown should not be used for purposes other than the production of permanent single crowns, inlays, onlays, and veneers. Any deviation from these Instructions For Use can have negative effects on the chemical and physical quality of objects made from Permanent Crown.

### CAUTION

Federal law restricts this device to sale or on the order of a dentist. For professional use only - Rx Only.

## 3. Specific Manufacturing Considerations

### NOTIFICATION

The device specifications have been validated using the printer process parameters indicated below.

### REQUIREMENTS

All accessories used must be designated for Permanent Crown Resin only. For full compliance and biocompatibility, Permanent Crown Resin requires a dedicated resin tank, Stainless Steel Build Platform, Form Wash, and finishing kit which should not be mixed with any other resins.

## RECOMMENDED 3D PRINTER AND PRINTING PARAMETERS

- a. Hardware: Formlabs SLA 3D Printer - Form 3B, Form 3B+ Visit [dental.formlabs.com](https://dental.formlabs.com) for printer compatibility
- Laser wavelength: 405 nm
  - Form 3 Resin Tank. Visit [dental.formlabs.com](https://dental.formlabs.com) for resin tank compatibility
  - Stainless Steel Build Platform only
- b. Software: Formlabs PreForm
- STL file import
  - Manual/Automatic rotation and placement
  - Manual/Automatic generation of supports
- c. Printing Parameters:
- Layer thickness: 50 µm
  - Optimal orientation: horizontal orientation, occlusal plane facing the build platform
  - Crowns, inlays, onlays, veneers: Minimum wall thickness: 1.0 mm
- d. Recommended Post-Processing Equipment:
- Formlabs Form Wash
  - Isopropyl alcohol (IPA) ≥99%
  - Formlabs Form Cure (390-405 nm)
  - Sandblaster
  - Glass bead blasting material - Perlablast® micro\* ( BEGO Co.) recommended
  - Finishing kit
  - Dental composite polish and tools

Post curing steps are described in Section 5.A. Print and Post-processing.

## CAUTION

The device specifications have been validated using the software, printers, and process parameters specified in this document. Any unauthorized changes to the process equipment, parameters, or software may result in a device that is out of specification and not covered under the FDA clearance. Users shall follow this document in order to use the Permanent Crown Resin. Users shall also follow the Instructions For Use documents and all maintenance requirements for the equipment identified in this document. Contact Formlabs for a list of compatible components.

\*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

## 4. Hazards And Precautions

### HAZARDS

Permanent Crown Resin contains polymerizable monomers which may cause skin irritation, allergic skin reaction, serious eye irritation, respiratory irritation or other allergic reactions in susceptible persons. If resin contacts skin, wash thoroughly with soap and water.

If any discomfort or adverse health effects occur from contact with parts printed with Permanent Crown Resin, discontinue use of this material.

Detailed information about the handling of Permanent Crown Resin can be found in the Safety Data Sheet and [Formlabs.com](https://www.formlabs.com).

1. **Eye Contact:** Can cause serious eye irritation. If eye irritation persists, seek medical attention.
2. **Skin Contact:** May cause skin irritation or allergic skin reactions. If skin irritation or rash occurs, seek medical attention.
3. **Inhalation:** May cause respiratory irritation. Avoid breathing mist, vapors, and sprays of this resin.
4. **Ingestion:** Avoid ingesting uncured or cured resin. If accidentally ingested, call a POISON CENTER and seek medical attention.
5. **Protection:** Protective clothing, glasses, and nitrile gloves should be worn while handling Permanent Crown Resin.
6. **Environment:** May cause harmful effects to aquatic life. Avoid release into the environment.
7. **Disposal:** Dispose of contents and cartridge as per local regulations.

**Contains:** Esterification products of 4,4'-isopropylidiphenol, ethoxylated and 2-methylprop-2enoic acid, silanized dental glass, methyl benzoylformate, diphenyl(2,4,6-trimethylbenzoyl) phosphine oxide.

Total content of inorganic fillers (particle size 0.7 $\mu$ m) is 30-50% by mass.

### CAUTION

Wash the printed part with solvent in a properly ventilated environment with proper protective masks and gloves.

1. Expired or uncured Permanent Crown Resin shall be disposed in accordance with local regulations.
2. IPA shall be disposed in accordance with local regulations.

## 5. Procedure for fabricating dental restorations with Permanent Crown Resin

Shake well before using a new cartridge of Permanent Crown Resin. Permanent Crown Resin in the resin tank does not need to be mixed prior to printing, up to a month of non-use, as solids (fillers) in the material do not settle.

### A. PRINTING AND POST-PROCESSING

1. **Set up:** Insert resin cartridge, Stainless Steel Build Platform, and the compatible resin tank into the Formlabs 3D printer.
2. **Printing:**
  - a. Prepare a print job using PreForm software. Import a dental restoration STL file. Orient files horizontally with the occlusal plane facing the build platform and generate supports. For recommendations on print orientation and support placement, see the detailed application guide at [support.formlabs.com](https://support.formlabs.com).
  - b. Send print job to printer. Begin print by selecting print job from print menu. Follow any prompts or dialogs shown on the printer screen. Printer will automatically complete the print.
3. **Washing:** Restorations printed with the Permanent Crown Resin should be washed using the Form Wash as follows,
  - a. Remove the build platform with the printed parts from the printer and place it into a Form Wash with clean IPA ( $\geq 99\%$ ).
  - b. Set the Form Wash for 3 minutes. Make sure the parts are fully submerged in IPA when washing. Exceeding wash duration may affect dimensional accuracy and performance of printed parts over time.
  - c. Once parts are washed, inspect parts for any remaining uncured resin on print surfaces. Use an IPA-filled squeeze bottle to remove any remaining resin on the printed parts and in between supports and rafts. An IPA-soaked brush may be used to help with uncured resin removal.
4. **Drying:** Use compressed air to dry the parts. The dried parts will have a white, powdery coating on the print surfaces. If any wet, uncured resin is still present after drying, reuse the IPA-soaked brush and IPA-filled squeeze bottle to remove uncured resin and dry again.
5. **Part removal:** After drying, remove printed parts from the build platform by wedging the part removal tool under the print raft and rotating the tool. For detailed techniques, see the support material on [dental.formlabs.com](https://dental.formlabs.com).
6. **Post cure:** Restorations printed with Permanent Crown Resin can be post cured in the Form Cure. To maintain dimensional accuracy and biocompatibility, post curing must be done in two steps:

**Caution:** It is imperative to use the compatible post-curing device listed below and follow the recommended cycle. The post-curing process is a critical step to meet the biocompatibility and material performance specifications of the restoration. No substitutions or changes to the process are allowed.

- a. With the raft and supports still intact, place printed parts in the Form Cure with the raft side down.
- b. Form Cure parts at 60°C (140°F) for 20 minutes.
- c. Remove parts from Form Cure and use a handpiece with a cutting disc to separate supports and raft from printed parts.
- d. Carefully sandblast the printed part surfaces to remove the white, powdery coating using a sandblaster with glass bead blasting material 50 µm, such as Perlablast® micro\* (e.g. from BEGO Co. REF 46092/54302) sandblasting medium at a maximum blasting pressure of 1.5 bar.
- e. Check printed part fit and finish and contour parts as needed using a dental handpiece and carbide rotary burs.
- f. Place finished parts in the Form Cure and post cure again at 60°C (140°F) for an additional 20 minutes, with the occlusal plane facing upwards. For best results and full compliance, printed parts must be post cured twice as recommended. Any deviations from this protocol, may lead to unsatisfactory fit and affect biocompatibility.

**Caution:** If there is an interruption or failure in the post-curing device cycle, the printed object should be not used until it has cured under a full cycle. Check the post-curing device instructions for use for how to properly resolve the post-curing device condition and then repeat the post-curing cycle with the printed objects.

**Attention:** Inspect the dental restorations for cracks. Discard if any cracks or damage are detected.

## B. POLISHING

Post cured parts must be polished using pumice stone and polishing compound prior to use.

1. High gloss polishing can be achieved with commonly used dental composite polishers.
2. Avoid overheating the parts during polishing.
3. Optimal surface quality is achieved by polishing after the full post curing cycle.
4. If rough marks are left behind after support removal, sand and polish to smooth surfaces to improve patient comfort.

\*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

## C. CLEANING AND DISINFECTION

Fully cured restorations made with Permanent Crown Resin can be easily cleaned and disinfected with IPA or ethanol. Steam cleaning using e.g. Triton SLA\* (BEGO Co.), or MD 520\* (Dürr Dental Co.) impression disinfectants are also possible. Follow the manufacturer's instructions.

## D. LUTING AGENT COMPATIBILITY

The finished restoration can be attached to a patient's tooth using self-adhesive cements (e. g. RelyX Unicem\*, 3M Espe\*) or composite cement with a primer (e. g. Variolink Esthetic DC\* and Monobond Plus\*, Ivoclar Vivadent\*). Follow the Instructions For Use provided by the luting agent manufacturer.

## E. STORAGE

1. Store finished parts in closed, opaque or amber containers in a cool, dry place out of direct sunlight.
2. Store the resin cartridges at room temperature (approx. 22°C or 72°F). Do not exceed 28°C (82°F) or drop below 4°C (39°F) when in storage.
3. Avoid ignition sources.
4. Do not use resin cartridges past the expiration date on the cartridge label.

## F. DISPOSAL

Any cured resin is non-hazardous and may be disposed of in regular waste.

1. Follow clinic or facility protocols for waste that is considered as a biohazard.
2. Liquid resin should be disposed of in accordance with governmental regulations (community, regional, or national).
  - a. Contact a licensed professional waste disposal service to dispose of liquid resin.
  - b. As with all foreign substances, do not allow resin to enter in storm or sewer drainage systems.
  - c. Avoid release into the environment.
  - d. Contaminated packaging: Dispose of as unused product.

## 6. Material properties and scope of delivery

### PHYSICAL DATA

Density (liquid resin)	ca. 1.4 - 1.5 g/cm <sup>3</sup>
Viscosity (liquid resin)	2,500 - 6,000 mPa*s
Flexural Strength	(Cured object) 116 MPa

\*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

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Wilh. Herbst GmbH & Co. KG**
**Manufactured for: Formlabs Inc.**